cystistat®

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cystistat®

STERILE SODIUM HYALURONATE SOLUTION

INDICATION: For temporary replacement of the glycosaminoglycan (GAG) layer in the bladder.

DESCRIPTION AND COMPOSITION: Cystistat® is a clear colourless solution presented in a 50 ml glass vial containing 40 mg sodium hyaluronate.

The glycosaminoglycan (GAG) layer on the luminal surface of the bladder wall is believed to provide a protective barrier against microorganisms, carcinogens, crystals and other agents present in the urine and has been identified as the primary defence mechanism in protecting the transitional epithelium from urinary irritants¹. Deficiencies in this GAG layer of the bladder epithelium may destroy its barrier function and allow the adherence of bacteria, microcrystals, proteins and ions, or the movement of ionic and nonionic solute residues (i.e. urea) across the epithelium². Cystistat® has been developed to temporarily replenish the deficient GAG layer on the bladder epithelium. The active substance is a highly purified sodium salt of hyaluronic acid.

DIRECTIONS: Instill the entire volume of this solution into the bladder after any residual urine has been removed. Instillation is intended to be performed by a qualified healthcare professional with the use of a sterile catheter and sterile syringe using appropriate hygienic technique, to minimise risk of infection. The sterile catheter and syringe are not supplied with the device. For best results Cystistat® should be retained in the bladder for as long as possible (a minimum of 30 minutes). Discard any unused portion.

There is evidence that the GAG layer of the bladder is deficient in cystitis. This deficiency contributes to the clinical symptoms² in the diseases such as interstitial cystitis³, cystitis caused by infections, trauma, urolithiasis, urinary retention, neoplasia and radiation induced cystitis. To alleviate cystitis associated with these conditions, it is recommended that Cystistat® be instilled into the bladder each week for 4-12 treatments and then monthly until symptoms resolve. The attending physician, urologist or radiologist should direct any prophylactic use of Cystistat®.

PRECAUTION: Do not administer to patients with known hypersensitivity reactions to sodium hyaluronate. Cystistat® is well tolerated, however, allergic reaction (local rash and/or itching) can occur. Catheterisation can also result in a temporary flare of symptoms, bladder/urethral pain and discomfort, or urinary tract infection. Inform your healthcare professional and discontinue use if any adverse reactions are experienced. As no clinical evidence is available on the use of sodium hyaluronate in children, pregnant and lactating women, treatment with Cystistat® is not recommended in these patients.

WARNING: Prior to administration, the user must examine the package to ensure the device is not damaged and is within the assigned expiry date. Inspect the solution for particulate matter and discoloration prior to administration.

DO NOT USE IF VIAL OR SEAL ARE DAMAGED.

DO NOT USE IF SOLUTION IS DISCOLOURED / CONTAINS PARTICULATE MATTER.

DO NOT USE AFTER THE EXPIRY DATE.

TO BE ADMINISTERED BY QUALIFIED HEALTHCARE PROFESSIONAL ONLY.

KEEP OUT OF THE REACH OF CHILDREN.

STORAGE: Store at room temperature (15-30°C). Do not freeze.

SUPPLIED: 1 x 50 ml vial of Cystistat®.

For single use only. Discard after use.

References:

- Messing, E.M. Interstitial cystitis and related syndromes. In: *Campbell's Urology*, 6th edition. Edited by P. C. Walsh, A. B. Retik, T. A. Stamey and E.D. Vaughan, Jr., Philadelphia: W. B. Saunders Co., Vol.1, Sect VII, Chapter 24, pp 982-1005, 1992.
- Parsons, C.L., D. Boychuk, S. Jones, R. Hurst and H. Callahan. Glycosaminoglycans: An Epithelial Permeability Barrier. *J.Urol.* 143: 139-142, 1990.
- Morales, A., Emerson, L., Nickel, J.C. and Lundie, M.: Intravesical hyaluronic acid in the treatment of refractory interstitial cystitis. *J. Urol* 156: 45-48, 1996

Further clinical information can be found at: www.cystistat.com

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© Mylan Teoranta, Ireland.